

Viruclear RC Series

Virus Removal Filter
Suitable for Complex Products



Viruclear RC Series

Virus Removal Filter Suitable for Complex Products



▲
Viruclear RC Silicone Cassettes
(Suspended Screen)

Virus filtration is a robust virus clearance method in the production process of common biopharmaceutical products. Due to its simple and gentle operation, clear mechanism, and easy validation, it has been widely used and can greatly improve the viral safety of biological products. The interception pore size of the Viruclear RC series is nominally 20nm, which is based on a size exclusion mechanism and can effectively remove various virus including parvovirus (the smallest virus currently known in nature).

In order to deal with products that are difficult to filter, Cobetter has launched a new series of virus removal filter made of Regenerated Cellulose. The natural and excellent hydrophilicity of regenerated cellulose materials enables lower adsorption, higher protein recovery and better filtration performance in complex products containing hydrophobic proteins, polymers and many impurities. The Viruclear RC series will be used as another solution to deal with the poor performance of the Viruclear VF Plus series. It is especially suitable for virus removal filtration in the extraction and preparation of biological products from human plasma, providing customers with a more comprehensive solution.

Features

- Unique surface hydrophilic improvement of PES membrane , which gives high mass capacity and high yield of virus retentive membrane
- Unique surface hydrophilic improvement of PES membrane, which gives high mass capacity and high yield of virus retentive membrane
- Composite supported membrane structure effectively improves the strength of the filter membrane
- Narrow pore size distribution gives the membrane good separation performance
- Robust virus removal capabilities
- Chemically stable
- Easy to install, use and test

Typical Applications

- Blood clotting factor
- Immunoglobulin
- Biochemical extracts
- Fusion protein
- Antibody fragments

Viruclear RC Series


Regenerated Cellulose Material Virus Removal Filter

Robust Removal of Virus

Test Condition

The test temperature is 21°C, set up six groups of parallel experiments for PP7 phage retention test to evaluate the parvovirus clearance ability.



| | |
|-----------------------|--|
| Filter | Viruclear RC |
| Model virus | PP7, added amount $>10^7$ PFU/mL |
| Load capacity up to | 50L/m ² , 100L/m ² , 150L/m ² , 200L/m ² |
| Pressure interruption | 30min |
| Post flush | 20L/m ² |



The standard solution selected for the test is 0.1g/L 0.22µm pre-filtered IVIG solution. Model virus selection: Pseudomonas aeruginosa PP7 phage is designated as a parvovirus model virus in the PDA TR41 rules. Viruclear RC products provide powerful virus removal capabilities and can handle the harshest conditions in the production process - multiple pressure interruptions. Viruclear RC maintains strong virus retention after multiple interruptions.

Linear Scalability

There is almost no significant difference in the water flux of Viruclear RC of different specifications (25°C, 29psi), and it maintains robust virus removal capabilities and has good linear amplification.




Viruclear RC H Series

RC material Virus Removal Filter Suitable for High Concentration Products

Viruclear RC H is a newly launched virus removal filter based on RC composite membrane. It further improves membrane strength and structural stability. It can maintain good pore size and flow channel structure, and can maintain high process flux when filtering high concentration protein solution, further improving the filtration ability for complex products and high concentration products while maintaining high efficiency virus removal ability. It is especially suitable for virus filtration of high concentration IVIG solution etc. in blood products industry.

Higher Process Throughput

A comparison of the filtration process flux of the Viruclear RC and Viruclear RC H filters, performed at room temperature and 30 psi under constant pressure filtration conditions, with the same 50 g/L intravenous gammaglobulin (IVIG), illustrates the difference in performance between the two filters when filtering high concentration proteins.



Comparison of flux decay trends over time for Viruclear RC and Viruclear RC H in series with Viruclear PNY for filtration of the same solution.

Robust Removal of Virus


Test Condition

Multiple batches of syringe filters are used for PP7 phage challenge experiments at room temperature, 30 psi.

| | |
|-------------------|---|
| Filter | Viruclear RC H Syringe Filter |
| Model virus | PP7, added amount >10 ⁷ PFU/mL |
| Feed solution | 0.1 g/L IVIG, Loading volume: 400L/m ² , Buffer flush volume: 20L/m ² |
| Constant pressure | 30 psi |




Test results




Virus Removal Filtration Process Development Service

Our technical engineers will work with customers to optimize the process parameters during the process development of virus removal filtration to obtain a robust, efficient and economical filtration operation.

Virus Removal Filterability Testing Device



Virus Removal Assembly Sketch Map



Optimization Strategy of Virus Removal Filtration Process

- Replace the traditional microporous prefilter with the adsorption-type nylon membrane filter PNY.
- Adopt Viruclear PDT depth prefilter.
- Increase the effective filtration area of PNY or PDT.
- Dilute, and reduce protein concentration.
- Optimize pH and conductivity of solution, and add stabilizer.

Product Specification

| | Specification | Effective Filtration Area | Processing Capacity |
|---|------------------|--|--|
| Viruclear RC Virus Removal Filter Series | Syringe filter | 2.5cm ² | Process development and virus removal validation |
| | Pilot scale | 0.015 / 0.08 / 0.3 / 0.5m ² | ≤ 100L |
| | Production scale | 1.50m ² | > 100L |

Product Specification

| | Specification | Membrane | Format | Shell | Effective Filtration Area | Assessries material |
|----------------------|------------------|-----------------------------------|--|-------------|---|---------------------|
| Virus Removal Filter | DS | RC | Syringe filter | PP | 2.5cm ² | Silicone |
| | Pilot scale | RC | Silicone cassette | Silicone | 0.015/0.08/ 0.3/0.5m ² | PP |
| | Production scale | RC | Silicone cassette | Silicone | 1.50m ² | PP |
| PNY Nylon Prefilter | DS series | Nylon | Syringe filter | PP | 3.4cm ² | Silicone |
| | Pilot scale | Nylon | Single-use Plastic-housing Cassette | PP | 0.025/0.12/0.30m ² | Silicone |
| | Production scale | Nylon | Single-use Plastic-housing Cassette | PP | 0.60/1.20/1.80m ² | Silicone |
| PDT Depth Prefilter | DS series | Cellulose, Diatomite, Nylon | Syringe filter | PP | 4.5cm ² | PP |
| | Pilot scale | Cellulose, Diatomite, Nylon | Single-use Plastic-housing Cassette | PP | 0.027/0.15/0.4m ² | PP |
| | Production scale | Cellulose, Diatomite, Nylon | Single-use Plastic-housing Cassette Single-use Plastic-housing Cassette (M02) | PP PP,GF | 0.92/1.1m ² 0.11/0.55/1.1m ² | PP PP |

Ordering Information

Viruclear PDT Virus Removal Prefilter (Depth Filter)

V P

Application

VP Viruclear PDF



D T

Material of Membrane

DT Depth filter

Type

DS Syringe filter
(4.5cm²)

Quantity/Package

N9 Only for syringe filter,
9pcs / pk
N1 1pcs / pk

P

Market

P Biopharmaceutical

Viruclear PDT Virus Removal Prefilter (Depth Filter)

V P

Application

VP Viruclear PDF



D T

Material of Membrane

DT Depth filter

Type

L4 L08TT(0.027m²)
SA CSCD(0.15m²)
SC CSCE(0.4m²)

Quantity/Package

N1 1pcs / pk

P

Market

P Biopharmaceutical

Viruclear PDT Virus Removal Prefilter (Depth Filter)

V P

Application

VP Viruclear PDF



D T

Material of Membrane

DT Depth filter

Type

SB CSCB(0.92m²)
SM CSCM(1.1m²)
M1 M02 Cassette(0.11m²)
M5 M02 Cassette(0.55m²)
MX M02 Cassette(1.1m²)

Quantity/Package

N1 1pcs / pk

P

Market

P Biopharmaceutical

Ordering Information

Viruclear PNY Virus Removal Prefilter (Nylon Filter)

V P

Application

VP Viruclear PNY

N Y

Material of Membrane

NY mNylon,
Modified nylon
membrane

Type

DS Syringe filter
(3.4cm²)

Quantity/Package

N9 Only for syringe filter,
9pcs / pk
N1 1pcs / pk

P

Market

P Biopharmaceutical



Viruclear PNY Virus Removal Prefilter (Nylon Filter)

V P

Application

VP Viruclear PNY

N Y

Material of Membrane

NY mNylon,
Modified nylon
membrane

Type

C02 250cm²
L02 0.12m²
L05 0.30m²

Type

TT 3/4" TC

Quantity/Package

N1 1pcs / pk

P

Market

P Biopharmaceutical



Viruclear PNY Virus Removal Prefilter (Nylon Filter)

V P

Application

VP Viruclear PNY

N Y

Material of Membrane

NY mNylon,
Modified nylon
membrane

Type

L10 0.60m²
L20 1.20m²
L30 1.80m²

Type

SS 1 1/2" TC

Quantity/Package

N1 1pcs / pk

P

Market

P Biopharmaceutical



Ordering Information

Viruclear RC Removal Syringe Filter with PP Shell (DS)

V F

Application

VF Virus filtration



R C

Material of Membrane

RC

Type

DS Syringe filter
(2.5cm²)

Quantity/Package

N9 Only for syringe filter,
9pcs / pk
(3 batches, 3 pcs per batch)
N1 1pcs / pk

P

Market

P Biopharmaceutical

Viruclear RC Silicone Cassette (Suspended Screen)

V F

Application

VF Virus filtration



R C

Material of Membrane

RC

Single-use Deflector

D Integration
N Non-integrated

Type

A1 Pilot module 1
Filtration Area 0.015m²
A2 Pilot module 2
Filtration Area 0.08m²
B1 Production module 1
Filtration Area 0.3m²
B2 Production module 2
Filtration Area 0.5m²
B3 Production module 3
Filtration Area 1.5m²

Quantity/Package

N1 1pcs / pk

P

Market

P Biopharmaceutical

Viruclear RC H Removal Syringe Filter with PP Shell (DS)

V F

VF Virus filtration



H

R C

Material of Membrane

RC

Type

DS Syringe filter
(2.5cm²)

Quantity/Package

N9 9 pcs / pk
N1 1pcs / pk

P

Market

P Biopharmaceutical

Viruclear RC H Silicone Cassette (Suspended Screen)

V F

Application

VF Virus filtration



H

R C

Material of Membrane

RC

Type

A1 Pilot module 1
Filtration Area 0.015m²
A2 Pilot module 2
Filtration Area 0.08m²
B1 Production module 1
Filtration Area 0.3m²
B2 Production module 2
Filtration Area 0.5m²
B3 Production module 3
Filtration Area 1.5m²

Quantity/Package

N1 1pcs / pk

P

Market

P Biopharmaceutical



Please contact us for more information
Hangzhou Cobetter Filtration Equipment Co.,Ltd.

Sales Add 13-18 F, Building 1, Zicheng International Innovation Center,
No.39 Jincheng Road, Xiaoshan District, Hangzhou 311215, China
Factory Cobetter Park, Heshang New Material Industrial Park,
Xiaoshan District, Hangzhou 311265, China

Tel +86 400-070-4266
Fax +86-571-87704256
www.cobetterfiltration.com